



30 Biowick fast response biohazard spill kits

BIOWICK

Dripwick | Biohazard Spill Kit | Maintenance

Base Fibre: granular



Product Code	Format	Size	Quantity in box
BIOWICK	Maintenance	250 gms	30/Box

Biowick is a fast response biological hazard spill kit tailored to deal with unwelcome spills such as blood, urine or vomit. The smartly designed box instantly converts into a scoop and receptacle. Bus lines, trains and schools are particularly suitable for the product as it can also be stored almost anywhere. There is nothing on the market that compares to this product.

Notes

- Ensure that disposal of these used sorbents complies with all regulations.

Spill kit's contents include:

250 grams of wick powder

Simply sprinkle the powder over the spill then collect it using the two halves of the scoop which converts from the containing box.

1 x pair of gloves

Waterproof gloves help to ensure safe handling of spilled substance.

2 x clean up wipes

Two clean up wipes are included for deodorising and finishing the final cleaning up process to leave the surface perfectly clean and dry.

1 x biohazard bag

The biohazard bag is used to contain and dispose of all the used items within the spill kit.



Biowick Fast Response spill kits
See this product in action on our website www.esp-us.com/Demo-Videos

* Basis weights and recycled content are ± 5%
 ** Measures a case/bale using ASTM 726 (long test) with medium viscosity oil and tap water. Statistical average derived from real time production QC data.
 *** Weight per case ± 0.5 kg/lbs.

Compliance

Related Regulations	Description
29CFR 1910.22(a)(2)	"The floor of every workroom shall be maintained in a clean and, so far as possible, a dry condition."
29CFR 1910.120(j)(1)(vii)	"... suitable quantities of proper absorbents shall be kept available and used in areas where spills, leaks or ruptures may occur."